



TrueKast

Standardized Offloading and Protection Systems

(patented technology)

Total Foot Protection is included in all Systems



- **STANDARDIZATION** of Cushioning, Splinting, and Casting Systems
- **CONVIENIENCE** of having each system in a single use kit.
- **KNOWLEDGE** of understanding what the patient will receive when ordered by the clinician.
- **CONSISTENCY** of Cushioning, Splinting, and Casting Methods.

TrueKast Total Cushioning System



Are you still using standard cast padding, gauze, and elastic wraps to protect wounds on the foot & ankle? Protect wounds on the foot & ankle, expensive skin substitutes or grafts, or the fragile skin of a newly closed wound while it regains its strength. Utilize the TCS for Diabetic Foot Ulcers that have healed while the patients wait for their custom orthotic or specialty shoe.

TrueKast Splinted Cushioning System



Advanced Cushioning Technology featuring “Reactive” Materials. Patent-Pending, Foam layer that lasts longer and works better than ordinary off-the-shelf dressings.

All-in-one kit designed for a Standardized approach to Offloading, Limited Fixation and Custom 7-Layer Fiberglass Posterior Splint.


TrueKast Total Contact Casting System



The TrueKast TCC system has a clinician-driven design providing one of the lightest, strongest, and most versatile TCC available today.

With an application time of less than eight minutes and a weight-bearing dry time of 8 minutes, the TrueKast provides the “Gold Standard” in offloading diabetic foot ulcers.

TrueKast QuickSaw Technology



The Integrated QuickSaw Cast Removal System provides for safe, fast, and patient-friendly removal of the Total Contact Cast while reducing the risks of patient injury due to the use of a powered, oscillating cast saw. It also reduces the amount of cast debris and aerosolized fiberglass cast dust produced by a vibrating power tool.

Everything Begins with Total Foot Protection